

## Recombinant Human TUBA1C

Catalog No: #GP12170



Package Size: #GP12170-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TUBA1C
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human TUBA1C
Target Name	tubulin, alpha 1c
Other Names	TUBA6; bcm948
Accession No.	Swissprot:Q9BQE3Gene Accession:BC004949
Uniprot	Q9BQE3
GeneID	84790;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only