

Recombinant Human TUSC2

Catalog No: #GP12172



Package Size: #GP12172-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TUSC2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	tumor suppressor candidate 2
Other Names	PAP; FUS1; PDAP2; C3orf11
Accession No.	Swissprot:O75896Gene Accession:BC023976
Uniprot	O75896
GeneID	11334;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only