

Recombinant Human TWIST2

Catalog No: #GP12173



Package Size: #GP12173-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TWIST2
Brief Description	Recombinant Protein
Immunogen Description	Full length of human TWIST2(172AA)
Target Name	twist family bHLH transcription factor 2
Other Names	FFDD3; DERMO1; SETLSS; bHLHa39
Accession No.	Swissprot:Q8WVJ9Gene Accession:BC017907
Uniprot	Q8WVJ9
GeneID	117581;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a basic helix-loop-helix type transcription factor and shares similarity with Twist. This protein may inhibit osteoblast maturation and maintain cells in a preosteoblast phenotype during osteoblast development. This gene may be upregulated in certain cancers. Mutations in this gene cause focal facial dermal dysplasia 3, Settleis type. Two transcript variants encoding the same protein have been found.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only