

Recombinant Human UACA

Catalog No: #GP12178



Package Size: #GP12178-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UACA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human uveal autoantigen with coiled-coil domains and ankyrin repeats
Target Name	uveal autoantigen with coiled-coil domains and ankyrin repeats
Other Names	NUCLING
Accession No.	Swissprot:Q9BZF9Gene Accession:BC113407
Uniprot	Q9BZF9
GenID	55075;
Storage	-20~-80°C, pH 7.6 PBS

Background

Regulates APAF1 expression and plays an important role in the regulation of stress-induced apoptosis. Promotes apoptosis by regulating three pathways, apoptosome up-regulation, LGALS3/galectin-3 down-regulation and NF-kappa-B inactivation. Regulates the redistribution of APAF1 into the nucleus after proapoptotic stress. Down-regulates the expression of LGALS3 by inhibiting NFKB1 (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only