

# Recombinant Human UBAP1

Catalog No: #GP12179



Package Size: #GP12179-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human UBAP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human UBAP1
Target Name	ubiquitin associated protein 1
Other Names	UAP; UBAP; NAG20; UBAP-1
Accession No.	Swissprot:Q9NZ09Gene Accession:BC098141
Uniprot	Q9NZ09
GeneID	51271;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is a member of the UBA domain family, whose members include proteins having connections to ubiquitin and the ubiquitination pathway. The ubiquitin associated domain is thought to be a non-covalent ubiquitin binding domain consisting of a compact three helix bundle. This particular protein originates from a gene locus in a refined region on chromosome 9 undergoing loss of heterozygosity in nasopharyngeal carcinoma (NPC).?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only