

Recombinant Human UBL4A

Catalog No: #GP12183



Package Size: #GP12183-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UBL4A
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ubiquitin-like 4A
Other Names	GDX; G6PD; GET5; MDY2; UBL4; TMA24; DX254E; DXS254E
Accession No.	Swissprot:P11441Gene Accession:BC043346
Uniprot	P11441
GeneID	8266;
Storage	-20~-80°C, pH 7.6 PBS

Background

UBL4A, also known as GDX, has been characterized as a ubiquitously expressed and highly conserved protein from humans to plants. The N-terminal region of UBL4A shares 43% sequence identity with ubiquitin, however the C-terminal region has no homology to ubiquitin and it is therefore unlikely that UBL4A plays any role in targeting cellular proteins for degradation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only