

Recombinant Human UCHL3

Catalog No: #GP12186



Package Size: #GP12186-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UCHL3
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ubiquitin carboxyl-terminal esterase L3 (ubiquitin thiolesterase)
Other Names	UCH-L3
Accession No.	Swissprot:P15374Gene Accession:BC018125
Uniprot	P15374
GeneID	7347;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the deubiquitinating enzyme family. Members of this family are proteases that catalyze the removal of ubiquitin from polypeptides and are divided into five classes, depending on the mechanism of catalysis. This protein may hydrolyze the ubiquitinyl-N-epsilon amide bond of ubiquitinated proteins to regenerate ubiquitin for another catalytic cycle. Alternative splicing results in multiple transcript variants that encode different protein isoforms.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only