

Recombinant Human UGGT2

Catalog No: #GP12187



Package Size: #GP12187-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UGGT2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	UDP-glucose glycoprotein glucosyltransferase 2
Other Names	UGT2; HUGT2; UGCGL2
Accession No.	Swissprot:Q9NYU1 Gene Accession:BC032302
Uniprot	Q9NYU1
GeneID	55757;
Storage	-20~-80°C, pH 7.6 PBS

Background

UDP-glucose:glycoprotein glucosyltransferase (UGT) is a soluble protein of the endoplasmic reticulum (ER) that selectively reglucosylates unfolded glycoproteins, thus providing quality control for protein transport out of the ER.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only