

# Recombinant Human PPP5C

Catalog No: #GP12192



Package Size: #GP12192-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PPP5C
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human PPP5C
Target Name	protein phosphatase 5, catalytic subunit
Other Names	PP5; PPT; PPP5
Accession No.	Swissprot:P53041?Gene Accession:BC001970
Uniprot	P53041
GeneID	5536;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a serine/threonine phosphatase which is a member of the protein phosphatase catalytic subunit family. Proteins in this family participate in pathways regulated by reversible phosphorylation at serine and threonine residues; many of these pathways are involved in the regulation of cell growth and differentiation. The product of this gene has been shown to participate in signaling pathways in response to hormones or cellular stress, and elevated levels of this protein may be associated with breast cancer development.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only