

Recombinant Human ALDH3A2

Catalog No: #GP12195



Package Size: #GP12195-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ALDH3A2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human ALDH3A2
Target Name	aldehyde dehydrogenase 3 family, member A2
Other Names	SLS; FALDH; ALDH10
Accession No.	Swissprot:P51648Gene Accession:BC002430
Uniprot	P51648
GeneID	224;
Storage	-20~-80°C, pH 7.6 PBS

Background

Aldehyde dehydrogenase isozymes are thought to play a major role in the detoxification of aldehydes generated by alcohol metabolism and lipid peroxidation. This gene product catalyzes the oxidation of long-chain aliphatic aldehydes to fatty acid. Mutations in the gene cause Sjogren-Larsson syndrome. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only