

# Recombinant Human SWAP70

Catalog No: #GP12198



Package Size: #GP12198-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SWAP70
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human SWAP70
Target Name	SWAP switching B-cell complex 70kDa subunit
Other Names	HSPC321; SWAP-70
Accession No.	Swissprot:Q9UH65Gene Accession:BC000616
Uniprot	Q9UH65
GeneID	23075;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Phosphatidylinositol 3,4,5-trisphosphate-dependent guanine nucleotide exchange factor (GEF) which, independently of RAS, transduces signals from tyrosine kinase receptors to RAC. It also mediates signaling of membrane ruffling. Regulates the actin cytoskeleton as an effector or adapter protein in response to agonist stimulated phosphatidylinositol (3,4)-bisphosphate production and cell protrusion (By similarity).

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only