

Recombinant Human UXT

Catalog No: #GP12211



Package Size: #GP12211-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UXT
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ubiquitously-expressed, prefoldin-like chaperone
Other Names	STAP1; ART-27
Accession No.	Swissprot:Q9UBK9Gene Accession:BC000720
Uniprot	Q9UBK9
GeneID	8409;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene functions as a cofactor that modulates androgen receptor-dependent transcription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this gene may play a role in tumorigenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only