

Recombinant Human VAMP4

Catalog No: #GP12212



Package Size: #GP12212-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VAMP4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	vesicle-associated membrane protein 4
Other Names	VAMP-4; VAMP24
Accession No.	Swissprot:O75379Gene Accession:BC007019
Uniprot	O75379
GeneID	8674;
Storage	-20~-80°C, pH 7.6 PBS

Background

Synaptobrevins/VAMPs, syntaxins, and the 25-kD synaptosomal-associated protein SNAP25 are the main components of a protein complex involved in the docking and/or fusion of synaptic vesicles with the presynaptic membrane. The protein encoded by this gene is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family. This protein may play a role in trans-Golgi network-to-endosome transport.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only