

Recombinant Human VAMP5

Catalog No: #GP12213



Package Size: #GP12213-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VAMP5
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	vesicle-associated membrane protein 5
Accession No.	Swissprot:O95183Gene Accession:BC017891
Uniprot	O95183
GeneID	10791;
Storage	-20~-80°C, pH 7.6 PBS

Background

Synaptobrevins/VAMPs, syntaxins, and the 25-kD synaptosomal-associated protein are the main components of a protein complex involved in the docking and/or fusion of vesicles and cell membranes. The VAMP5 gene is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family and the SNARE superfamily. This VAMP family member may participate in vesicle trafficking events that are associated with myogenesis.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only