

# Recombinant Human VAPA

Catalog No: #GP12215



Package Size: #GP12215-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human VAPA
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa
Other Names	VAP-A; VAP33; VAP-33; hVAP-33
Accession No.	Swissprot:Q9P0L0Gene Accession:BC002992
Uniprot	Q9P0L0
GeneID	9218;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a type IV membrane protein. It is present in the plasma membrane and intracellular vesicles. It may also be associated with the cytoskeleton. This protein may function in vesicle trafficking, membrane fusion, protein complex assembly and cell motility. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only