

Recombinant Human VASH1

Catalog No: #GP12216



Package Size: #GP12216-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VASH1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	vasohibin 1
Other Names	KIAA1036
Accession No.	Swissprot:Q7L8A9Gene Accession:BC009031
Uniprot	Q7L8A9
GeneID	22846;
Storage	-20~-80°C, pH 7.6 PBS

Background

Vasohibin-1 is a protein that in humans is encoded by the VASH1 gene. Angiogenesis inhibitor. Inhibits migration, proliferation and network formation by endothelial cells as well as angiogenesis. This inhibitory effect is selective to endothelial cells as it does not affect the migration of smooth muscle cells or fibroblasts. Does not affect the proliferation of cancer cells in vitro, but inhibits tumor growth and tumor angiogenesis. Acts in an autocrine manner. Inhibits artery neointimal formation and macrophage infiltration.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only