

# Recombinant Human VASH2

Catalog No: #GP12217



Package Size: #GP12217-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human VASH2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	vasohibin 2
Accession No.	Swissprot:Q86V25Gene Accession:BC028194
Uniprot	Q86V25
GeneID	79805;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Vasohibin-2, also known as VASH2, is a 355 amino acid protein belonging to the Vasohibin family. Expressed in various embryonic organs, Vasohibin-2 can be detected during embryonic weeks six through twelve and is found in vessels of 20-week embryonic organs as well as in endothelial cells of neonatal large vessels. Induced by VEGF, Vasohibin-2 is an angiogenesis inhibitor and interferes with the proliferation, migration and network formation by endothelial cells.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only