

Recombinant Human VCAM1

Catalog No: #GP12218



Package Size: #GP12218-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VCAM1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human VCAM1
Target Name	vascular cell adhesion molecule 1
Other Names	CD106; INCAM-100
Accession No.	Swissprot:P19320Gene Accession:BC068490
Uniprot	P19320
GeneID	7412;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the Ig superfamily and encodes a cell surface sialoglycoprotein expressed by cytokine-activated endothelium. This type I membrane protein mediates leukocyte-endothelial cell adhesion and signal transduction, and may play a role in the development of atherosclerosis and rheumatoid arthritis. Three alternatively spliced transcripts encoding different isoforms have been described for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only