

Recombinant Human VIL1

Catalog No: #GP12221



Package Size: #GP12221-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VIL1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	villin 1
Other Names	VIL; D2S1471
Accession No.	Swissprot:P09327Gene Accession:BC017303
Uniprot	P09327
GeneID	7429;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of a family of calcium-regulated actin-binding proteins. This protein represents a dominant part of the brush border cytoskeleton which functions in the capping, severing, and bundling of actin filaments. Two mRNAs of 2.7 kb and 3.5 kb have been observed; they result from utilization of alternate poly-adenylation signals present in the terminal exon.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only