

## Recombinant Human VPS25

Catalog No: #GP12227

Package Size: #GP12227-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human VPS25
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	vacuolar protein sorting 25 homolog (S. cerevisiae)
Other Names	DERP9; EAP20; FAP20
Accession No.	Swissprot:Q9BRG1Gene Accession:BC006282
Uniprot	Q9BRG1
GeneID	84313;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a protein that is a subunit of the endosomal sorting complex required for transport II (ESCRT-II). This protein complex functions in sorting of ubiquitinated membrane proteins during endocytosis. A pseudogene of this gene is present on chromosome 1.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only