

Recombinant Human WDR61

Catalog No: #GP12244



Package Size: #GP12244-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human WDR61
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	WD repeat domain 61
Other Names	SKI8; REC14
Accession No.	Swissprot:Q9GZS3Gene Accession:BC010080
Uniprot	Q9GZS3
GeneID	80349;
Storage	-20~-80°C, pH 7.6 PBS

Background

WDR61 is a subunit of the human PAF and SKI complexes, which function in transcriptional regulation and are involved in events downstream of RNA synthesis, such as RNA surveillance (Zhu et al., 2005 [PubMed 16024656]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only