

## Recombinant Human XAB2

Catalog No: #GP12256



Package Size: #GP12256-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human XAB2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human XAB2
Target Name	XPA binding protein 2
Other Names	HCNP; HCRN; SYF1; NTC90
Accession No.	Swissprot:Q9HCS7Gene Accession:BC007208
Uniprot	Q9HCS7
GeneID	56949;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Involved in transcription-coupled repair (TCR), transcription and pre-mRNA splicing.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only