

Recombinant Human YY1

Catalog No: #GP12259



Package Size: #GP12259-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human YY1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 215-414 amino acids of human YY1
Target Name	YY1 transcription factor
Other Names	DELTA; NF-E1; UCRBP; INO80S; YIN-YANG-1
Accession No.	Swissprot:P25490Gene Accession:BC037308
Uniprot	P25490
GeneID	7528;
Storage	-20~-80°C, pH 7.6 PBS

Background

YY1 is a ubiquitously distributed transcription factor belonging to the GLI-Kruppel class of zinc finger proteins. The protein is involved in repressing and activating a diverse number of promoters. YY1 may direct histone deacetylases and histone acetyltransferases to a promoter in order to activate or repress the promoter, thus implicating histone modification in the function of YY1.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only