

Recombinant Human ZFAND3

Catalog No: #GP12264



Package Size: #GP12264-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZFAND3
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	zinc finger, AN1-type domain 3
Other Names	TEX27
Accession No.	Swissprot:Q9H8U3Gene Accession:BC031481
Uniprot	Q9H8U3
GeneID	60685;
Storage	-20~-80°C, pH 7.6 PBS

Background

ZFAND3 contains 1 A20-type zinc finger and 1 AN1-type zinc finger. The exact function of ZFAND3 remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only