

# Recombinant Human ZFP28

Catalog No: #GP12266



Package Size: #GP12266-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ZFP28
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human ZFP28
Target Name	ZFP28 zinc finger protein
Other Names	mkr5
Accession No.	Swissprot:Q8NHY6Gene Accession:BC028416
Uniprot	Q8NHY6
GeneID	140612;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May be involved in transcriptional regulation. May have a role in embryonic development.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only