

# Recombinant Human ZFP2

Catalog No: #GP12269



Package Size: #GP12269-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ZFP2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human ZFP2
Target Name	ZFP2 zinc finger protein
Other Names	ZNF751
Accession No.	Swissprot:Q6ZN57Gene Accession:BC132781
Uniprot	Q6ZN57
GeneID	80108;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Probable transcription factor involved in neuronal differentiation and/or phenotypic maintenance.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only