

Recombinant Human ZNF23

Catalog No: #GP12271



Package Size: #GP12271-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF23
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 50-149 amino acids of human ZNF23
Target Name	zinc finger protein 23
Other Names	KOX16; ZNF359; ZNF612; Zfp612
Accession No.	Swissprot:P17027Gene Accession:BC007256
Uniprot	P17027
GeneID	7571;
Storage	-20~-80°C, pH 7.6 PBS

Background

May be involved in transcriptional regulation. May have a role in embryonic development.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only