

Recombinant Human ZNF239

Catalog No: #GP12272



Package Size: #GP12272-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF239
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human ZNF239
Target Name	zinc finger protein 239
Other Names	MOK2; HOK-2
Accession No.	Swissprot:Q16600Gene Accession:BC026030
Uniprot	Q16600
GeneID	8187;
Storage	-20~-80°C, pH 7.6 PBS

Background

MOK2 proteins are DNA- and RNA-binding proteins that are mainly associated with nuclear RNP components, including the nucleoli and extranucleolar structures (Arranz et al., 1997 [PubMed 9121460]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only