

Recombinant Human RNF115

Catalog No: #GP12274



Package Size: #GP12274-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RNF115
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ring finger protein 115
Other Names	BCA2; ZNF364
Accession No.	Swissprot:Q9Y4L5Gene Accession:BC054049
Uniprot	Q9Y4L5
GeneID	27246;
Storage	-20~-80°C, pH 7.6 PBS

Background

E3 ubiquitin-protein ligase that mediates E2-dependent, 'Lys-48'- and/or 'Lys-63'-linked polyubiquitination of substrates and may play a role in diverse biological processes. Through their polyubiquitination, may play a role in the endosomal trafficking and degradation of membrane receptors including EGFR, FLT3, MET and CXCR4.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only