

Recombinant Human TBCA

Catalog No: #GP12306



Package Size: #GP12306-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TBCA
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	tubulin folding cofactor A
Accession No.	Swissprot:O75347Gene Accession:BC018210
Uniprot	O75347
GeneID	6902;
Storage	-20~-80°C, pH 7.6 PBS

Background

The product of this gene is one of four proteins (cofactors A, D, E, and C) involved in the pathway leading to correctly folded beta-tubulin from folding intermediates. Cofactors A and D are believed to play a role in capturing and stabilizing beta-tubulin intermediates in a quasi-native confirmation. Cofactor E binds to the cofactor D/beta-tubulin complex; interaction with cofactor C then causes the release of beta-tubulin polypeptides that are committed to the native state.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only