

Recombinant Human CARTPT

Catalog No: #GP12313



Package Size: #GP12313-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CARTPT
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	CART prepropeptide
Other Names	CART
Accession No.	Swissprot:Q16568Gene Accession:BC029882
Uniprot	Q16568
GeneID	9607;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a secreted protein which is processed by prohormone/proprotein convertases to produce smaller, biologically active peptides. Expression of the transcript for this gene is regulated by certain drugs such as cocaine, and the encoded protein is thought to be involved in the regulation of appetite and stress. Mutations in this gene are associated with susceptibility to obesity.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only