

Recombinant Human ERMP1

Catalog No: #GP12325



Package Size: #GP12325-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ERMP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 705-904 amino acids of human ERMP1
Target Name	endoplasmic reticulum metallopeptidase 1
Other Names	FXNA; KIAA1815; bA207C16.3
Accession No.	Swissprot:Q7Z2K6Gene Accession:BC031630
Uniprot	Q7Z2K6
GeneID	79956;
Storage	-20~-80°C, pH 7.6 PBS

Background

Within the ovary, required for the organization of somatic cells and oocytes into discrete follicular structures.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only