

Recombinant Human CBWD1

Catalog No: #GP12329



Package Size: #GP12329-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CBWD1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 196-395 amino acids of human CBWD1
Target Name	COBW domain containing 1
Other Names	COBP
Accession No.	Swissprot:Q9BRT8Gene Accession:BC132759
Uniprot	Q9BRT8
GeneID	55871;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only