

Recombinant Human MYBPHL

Catalog No: #GP12330



Package Size: #GP12330-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MYBPHL
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 155-354 amino acids of human MYBPHL
Target Name	myosin binding protein H-like
Accession No.	Swissprot:A2RUH7Gene Accession:BC132902
Uniprot	A2RUH7
GeneID	343263;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein with two immunoglobulin superfamily domains and a fibronectin 3 domain. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only