

Recombinant Human RUFY4

Catalog No: #GP12340



Package Size: #GP12340-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RUFY4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 372-571 amino acids of human RUFY4
Target Name	RUN and FYVE domain containing 4
Other Names	ZFYVE31
Accession No.	Swissprot:Q6ZNE9Gene Accession:BC113712
Uniprot	Q6ZNE9
GeneID	285180;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only