

Recombinant Human CDHR4

Catalog No: #GP12344



Package Size: #GP12344-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CDHR4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 18-180 amino acids of human CDHR4
Target Name	cadherin-related family member 4
Other Names	CDH29; PRO34300
Accession No.	Swissprot:A6H8M9Gene Accession:BC132756
Uniprot	A6H8M9
GeneID	389118;
Storage	-20~-80°C, pH 7.6 PBS

Background

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only