

Recombinant Human SLC35E2

Catalog No: #GP12345



Package Size: #GP12345-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SLC35E2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 187-266 amino acids of human SLC35E2
Target Name	solute carrier family 35, member E2
Accession No.	Swissprot:P0CK97Gene Accession:BC112289
Uniprot	P0CK97
GeneID	9906;
Storage	-20~-80°C, pH 7.6 PBS

Background

SLC35E2 is a putative transporter.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only