

# Recombinant Human PRSS38

Catalog No: #GP12347



Package Size: #GP12347-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PRSS38
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	protease, serine, 38
Other Names	MPN2
Accession No.	Swissprot:A1L453Gene Accession:BC130400
Uniprot	A1L453
GeneID	339501;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The function of this protein remains unknown.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only