

# Recombinant Human ITLN2

Catalog No: #GP12350



Package Size: #GP12350-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ITLN2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 125-324 amino acids of human ITLN2
Target Name	intelectin 2
Other Names	HL2; HL-2
Accession No.	Swissprot:Q8WWU7Gene Accession:BC117225
Uniprot	Q8WWU7
GeneID	142683;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May play a role in the defense system against pathogens.Expressed only in the small intestine.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only