

Recombinant Human OIT3

Catalog No: #GP12356



Package Size: #GP12356-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human OIT3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 122-321 amino acids of human OIT3
Target Name	oncoprotein induced transcript 3
Other Names	LZP
Accession No.	Swissprot:Q8WWZ8Gene Accession:BC126440
Uniprot	Q8WWZ8
GeneID	170392;
Storage	-20~-80°C, pH 7.6 PBS

Background

May be involved in hepatocellular function and development.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only