

Recombinant Human RSPH1

Catalog No: #GP12357



Package Size: #GP12357-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RSPH1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 160-309 amino acids of human RSPH1
Target Name	radial spoke head 1 homolog (Chlamydomonas)
Other Names	CT79; TSA2; RSP44; TSGA2; RSPH10A
Accession No.	Swissprot:Q8WYR4 Gene Accession:BC113367
Uniprot	Q8WYR4
GeneID	89765;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a male meiotic metaphase chromosome-associated acidic protein. This gene is expressed in tissues with motile cilia or flagella, including the trachea, lungs, airway brushings, and testes. Mutations in this gene result in primary ciliary dyskinesia-24. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only