

Recombinant Human RAB41

Catalog No: #GP12360



Package Size: #GP12360-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RAB41
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	RAB41, member RAS oncogene family
Accession No.	Swissprot:Q5JT25Gene Accession:BC117239
Uniprot	Q5JT25
GeneID	347517;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a small GTP-binding protein that belongs to the largest family within the Ras superfamily. These proteins function as regulators of membrane trafficking. They cycle between inactive GDP-bound and activated GTP-bound states, which is controlled by GTP hydrolysis-activating proteins (GAPs). This family member can be activated by the GAP protein RN-Tre, and it is localized to the Golgi complex.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only