

Recombinant Human COMMD6

Catalog No: #GP12363



Package Size: #GP12363-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human COMMD6
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	COMM domain containing 6
Other Names	Acrg
Accession No.	Swissprot:Q7Z4G1 Gene Accession:BC117391
Uniprot	Q7Z4G1
GeneID	170622;
Storage	-20~-80°C, pH 7.6 PBS

Background

COMMD6 belongs to a family of NF-kappa-B (see RELA; MIM 164014)-inhibiting proteins characterized by the presence of a COMM domain (see COMMD1; MIM 607238) (de Bie et al., 2006 [PubMed 16573520])

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only