

Recombinant Human ZNRD1

Catalog No: #GP12364



Package Size: #GP12364-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNRD1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	zinc ribbon domain containing 1
Other Names	TEX6; ZR14; HTEX6; Rpa12; hZR14; HTEX-6; TCTEX6; tctex-6
Accession No.	Swissprot:Q9P1U0Gene Accession:BC050608
Uniprot	Q9P1U0
GeneID	30834;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a DNA-directed RNA polymerase I subunit. The encoded protein contains two potential zinc-binding motifs and may play a role in regulation of cell proliferation. The encoded protein may be involved in cancer and human immunodeficiency virus progression. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only