

Recombinant Human FAM101A

Catalog No: #GP12370



Package Size: #GP12370-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FAM101A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 82-216 amino acids of human FAM101A
Target Name	family with sequence similarity 101, member A
Other Names	CFM2
Accession No.	Swissprot:Q6ZTI6Gene Accession:BC141805
Uniprot	Q6ZTI6
GeneID	100533183;144347;
Storage	-20~-80°C, pH 7.6 PBS

Background

Involved in the regulation of the perinuclear actin network and nuclear shape through interaction with filamins. Plays an essential role in actin cytoskeleton formation in developing cartilaginous cells.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only