

# Recombinant Human NUDT21

Catalog No: #GP12374



Package Size: #GP12374-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human NUDT21
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	nudix (nucleoside diphosphate linked moiety X)-type motif 21
Other Names	CPSF5; CFIM25
Accession No.	Swissprot:O43809Gene Accession:BC001403
Uniprot	O43809
GeneID	11051;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is one subunit of a cleavage factor required for 3' RNA cleavage and polyadenylation processing. The interaction of the protein with the RNA is one of the earliest steps in the assembly of the 3' end processing complex and facilitates the recruitment of other processing factors. This gene encodes the 25kD subunit of the protein complex, which is composed of four polypeptides.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only