

# Recombinant Human CNPY2

Catalog No: #GP12395



Package Size: #GP12395-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CNPY2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	canopy FGF signaling regulator 2
Other Names	MSAP; TMEM4; ZSIG9; HP10390
Accession No.	Swissprot:Q9Y2B0Gene Accession:BC065015
Uniprot	Q9Y2B0
GeneID	10330;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Positive regulator of neurite outgrowth by stabilizing myosin regulatory light chain (MRLC). It prevents MIR-mediated MRLC ubiquitination and its subsequent proteasomal degradation.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only