

# Recombinant Human FAM96A

Catalog No: #GP12401



Package Size: #GP12401-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human FAM96A
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	family with sequence similarity 96, member A
Accession No.	Swissprot:Q9H5X1 Gene Accession:BC008865
Uniprot	Q9H5X1
GeneID	84191;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May play a role in chromosome segregation through establishment of sister chromatid cohesion.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only