

# Recombinant Human LSAMP

Catalog No: #GP12407



Package Size: #GP12407-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human LSAMP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 29-315 amino acids of human LSAMP
Target Name	limbic system-associated membrane protein
Other Names	LAMP; IGLON3
Accession No.	Swissprot:Q13449Gene Accession:BC033803
Uniprot	Q13449
GeneID	4045;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a neuronal surface glycoprotein found in cortical and subcortical regions of the limbic system. During development of the limbic system, this encoded protein is found on the surface of axonal membranes and growth cones, where it acts as a selective homophilic adhesion molecule, and guides the development of specific patterns of neuronal connections.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only