

Recombinant Human DYDC2

Catalog No: #GP12410



Package Size: #GP12410-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DYDC2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	DPY30 domain containing 2
Accession No.	Swissprot:Q96IM9Gene Accession:BC007374
Uniprot	Q96IM9
GeneID	84332;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of a family of proteins that contains a DPY30 domain. This gene locus overlaps with a closely related gene on the opposite strand. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only